

Ensure your equipment operates successfully with NIST-Traceable Calibration Service



NIST-Traceable Calibration

Calibration is an essential part of equipment maintenance. STERIS offers calibration procedures to your equipment's specific needs over its entire lifecycle.

Performed by specially trained, field-based technicians, all calibration tests are fully NIST-traceable and detailed in a complete report kept on file with STERIS for a minimum of three (3) years.



In addition to NIST-Traceable Calibration, visit our website to learn about more services offered by **STERIS** at sterislifesciences.com

Calibration and verification points are determined by STERIS, the manufacturer, or your SOPs. We will make necessary adjustments to 'out of tolerance' conditions we may find and include the results in your final report.

A full report is offered after every calibration. Your STERIS technician will take the time to review your report that will contain:

- 1.) A list of 'As Found' and 'Corrected' readings
- 2.) Tolerance page and signature page including a list of standards, procedures, temperatures, pressures, comments, and signatures from both you and your STERIS technician
- 3.) Calibration certificates for each instrument used in the calibration

Your equipment will also be tagged to show cGMP compliance.



NIST-TRACEABLE CALIBRATION IS OFFERED ON:

- Sterilizers
- Washers
- Generators
- Water Stills
- Low temperature sterilization and biodecontamination (VHP) equipment



STERIS CORPORATION
5960 Heisley Road
Mentor, OH 44060, USA
T +1 800 444 9009

www.STERISLifeSciences.com

STERIS Offices Worldwide

Asia Pacific	+65 (0) 68 41 7677
Canada	+1 800 444 9009
Finland	+358 9 25851
France	+33 (0) 2 38 84 85 40
Germany	+49 (0) 221 466120
Italy	+39 (0) 2 21303 424
Latin America	+1 800 444 9009
Spain	+34 (0) 916 585 920
United Kingdom	+44 (0) 1256 840400

STERIS Life Sciences

For over 100 years, STERIS has been known as a global leader, trusted partner and solutions provider in the field of sterilization and contamination control. Today, STERIS continues building on this heritage by advancing the Science of sterilization, cleaning and infection control while offering Solutions that meet our Customers' needs and high standards. STERIS is dedicated to helping you enhance the Life of your patients and the life of your equipment. From Formulated Chemistries, to Capital Equipment, and to Parts and Services, STERIS Life Sciences is your Solutions provider.

© 2018 STERIS Corporation
All rights reserved

No part of this publication may be reproduced in any material form without the permission of STERIS Corp.